

Labelling DNA with EdU for flow cytometry detection - fission yeast

1. Label (tk+ ENT+) strain¹ with (1-10 μ M) EdU - can use YES or EMM. Cells seem to only incorporate for one cycle and show checkpoint arrest (i.e. only one strand can be labelled)².
2. Fix cells (OD600 0.1-0.5) by spinning down and resuspending in 70% etoh. Transfer to 1.5 ml eppendorf. Cells can be stored at this stage (4°C up to a week).
3. Rehydrate by adding cells (ca. 0.2 ml, about 2 x 10⁶ cells) to 5 ml PBS and spin down. Repeat wash PBS x 1 (this extra wash may not be necessary). If necessary split cells at this stage e.g. for sytox green determination of DNA content.
4. Spin down cells and pour off PBS

5. Assemble Click it reaction components as follows:

- 1X ClickIT reaction buffer (diluted from x10 stock in kit - keep @ 4°C)
- CuSO₄ (100mM)
- Fluorescent azide (kept at -20°C)
- 1X Reaction buffer additive (dilute from x10 stock in kit by diluting with deionized water)

6. Prepare appropriate volume of Click it cocktail as follows:

Table 3. Click-iT™ reaction cocktails.

Reaction components	Number of reactions						
	1	2	5	10	15	30	50
1X Click-iT™ Reaction Buffer (prepared in step 1.5)	438 μ L	875 μ L	2.19 mL	4.38 mL	6.57 mL	13.2 mL	21.9 mL
CuSO ₄ (Component H)	10 μ L	20 μ L	50 μ L	100 μ L	150 μ L	300 μ L	500 μ L
Fluorescent dye azide (prepared in step 1.3)	2.5 μ L	5 μ L	12.5 μ L	25 μ L	37.5 μ L	75 μ L	125 μ L
Reaction Buffer Additive (prepared in step 4.6)	50 μ L	100 μ L	250 μ L	500 μ L	750 μ L	1.5 mL	2.5 mL
Total volume	500 μ L	1 mL	2.5 mL	5 mL	7.5 mL	15 mL	25 mL

Add 0.5 ml click it cocktail to cell pellet from step 4. Vortex. Incubate in dark at room temp for 30 mins. Go to step 7

ALTERNATIVELY

5. Make up buffer containing 50 mM Tris, pH 7.4, 150 mM NaCl, 2 mM CuSO₄, 10 μ M Alexa Fluor 488 azide (added from 2mM stock in DMSO; A10266, Invitrogen), 10mM sodium ascorbate (added last from 0.1 M stock)

¹ Strains that are just TK+ (ie no hENT) can be used but a much higher conc of EdU must be used (see Hua and Kearsy 2011)

² EdU labelling seems to cause DNA damage and activates the DNA damage checkpoint, making labelling experiments over more than one cell cycle problematic. The kinetics of S phase appear not to be altered by EdU incorporation.

Labelling DNA with EdU for flow cytometry detection - fission yeast

6. Add 0.5 ml of this mix to cell pellet from step 4 and resuspend. Incubate in dark at room temp for 30 mins.

7. Spin down cells, discard supernatant and resuspend in 5 ml PBS.

8. Repeat PBS wash by spinning cells down again, and resuspending in 1ml PBS

9. Finally resuspend in 1 ml PBS

10. Transfer to flow cytometry tube. Sonicate to break up doublets (use Bioruptor bath sonicator, 20s, low power. **Do not use** Soniprep tip sonicator).

11. Sample is ready for analysis.

12. Samples can be examined by fluorescence microscopy by concentrating and mounting in 1.2% agarose (containing DAPI if required).

Manual from invitrogen:

<http://probes.invitrogen.com/media/pis/mp35002.pdf>

REFERENCE

Hua, H., and S.E. Kearsey. Monitoring DNA replication in fission yeast by incorporation of 5-ethynyl-2'-deoxyuridine. 2011. *Nucleic Acids Res.* 39:e60.