# MCM chromatin-binding assay

## January 2003

Stephen Kearsey Dept Zoology South Parks Rd Oxford

stephen.kearsey@zoo.ox.ac.uk

## **Buffers:**

1. EMM sorb pH 7.0

KH phallate	3g	15 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15 mM
NH <sub>4</sub> Cl	5 g	90 mM
Sorbitol	218.4 g	(1.2M)

pH with NaOH to 7.0

water to 11

2. Zymolyase 20-T (ICN 320921 20000U/g)

20mg/ml in EMM sorb buffer - store 50 µl aliquots at -70°C

source of zymolyase is important when digesting -N cells

3. Extraction buffer – depends on protein being examined. For Mcm2-7 use buffer 1 or buffer For Mcm2-7 in meiosis experiments use buffer 2 (no magnesium).

# **Extraction buffers for CBA assay**

## Extraction buffer 1 (normal +Mg2+)

As per buffer 2 but omit ETDA and add MgAc to 2 mM Also, do not add spermine and spermidine

## **Extraction buffer 2 (normal)**

for 500 ml
3.02 g
36.4
1 ml of 0.5 M
7.35 g

After autoclaving add to 50 ml lots in falcon tube

0.5 mM spermidine HCl	25 µl 1 M
0.15 mM spermine HCl	7.5 µl 1M

store at -20°C

# Extraction buffer 3 (high pH, low salt)

	for 500 ml
20 mM HEPES-KOH pH 7.7	
0.4 M sorbitol	36.4
1 mM EDTA	1 ml of 0.5 M
50 mM KAc	2.45 g
Autoclave	

After autoclaving add to 50 ml lots in falcon tube

0.5 mM spermidine HCl	25 µl 1 M
0.15 mM spermine HCl	7.5 µl 1M

store at -20°C

# **Extraction buffer 4 (ELB salts)**

for 500 ml

for 500 ml

3.02 g

36.4

2.5 mM MgCl2 50 mM KCl 10 mM HEPES-KOH ph 7.7

## Extraction buffer 5 (normal, low salt, for DNA polymerase alpha)

20 mM Pipes-KOH pH 6.8 0.4 M sorbitol

to 50 ml aliquots add

0.1 ml 5 M KAc	(final conc 10 mM)
0.5 mM spermidine HCl	25 µl 1 M
0.15 mM spermine HCl	7.5 μl 1M
0.1 ml 0.5 M EDTA pH 8.0	(final conc 1 mM)

#### 1. Basic procedure

- 1. Grow gfp-tagged strain strain into log phase in YE or EMM (EMM is recommended) (typically 20 ml @ OD<sub>595</sub> =0.2-0.5).
- 2. Label and place all tubes and solutions on ice.
- 3. Transfer around 20ml to Falcon tube. (*Add 500µl Na Azide*(*HIGHLY TOXIC*!- gloves, labcoat) (10%) (for eg time course where it is important to stop the cells in DNA replication).
- 4. Spin tubes  $(3k, 3 \min, 4^{\circ}C)$ .
- 5. Resuspend in 1.8ml EMM sorb buffer and transfer to 2ml eppendorf. (can include 10 mM DTT in this step for cells that may be difficult to digest eg after -n2 starvation).
- 6. Spin tubes  $(4k, 1 \min, 4^{\circ}C)$ .
- 7. Resuspend cells in 450µl EMM sorb buffer containing DTT (10 mM, add 1/100 1M stock).
- 8. Add 50 µl of 20 mg/ml zymolyase 20-T (ICN 320921 20000U/g).

Digest at 32°C for 10-20 minutes. May need to use a much higher concentration of zymolyase or longer time of digestion, if cells have been nitrogen starved eg in meiosis release experiments.

- 9. Test for adequate digestion by adding SDS to aliquot cells should go dark under phase (aim for >95% phase dark).
- 10. Add 1.5 ml EMM sorb buffer (at 0°C).
- 11. Spin tubes (4k, 1 min, 4°C). (Spin1)
- 12. Resuspend in 2 ml EMM sorb buffer (at 0°C).
- 13. Spin tubes (4k, 1 min, 4°C). (Spin2)
- 14. Resuspend in 2 ml extraction buffer (at 0°C).\*\*
- 15. Spin tubes (4k, 1 min, 4°C). (Spin3)
- 16. Resuspend in 0.9 ml extraction buffer **containing protease inhibitor** (Roche 1 836 153 "complete mini)\*.
- 17. Split into 2 x 450µl.
- 18. To one  $450\mu$ l add 50 µl extraction buffer plus 10% Triton X-100(made in extraction buffer) = +Triton. Remaining 0.45 ml serves as minus Triton control).
- 19. Mix both tubes by inverting. Transfer +Triton to 20°C water bath. Keep -Triton on ice.
- 20. Incubate for 5 min, inverting periodically. Spin tubes (4k, 1 min, 4°C).
- 21. Take off supernatant, flick tubes so cells come off of tube wall and add 1.5 ml methanol.
- 22. Leave on ice for at least 6 mins.
- 23. Spin tubes (4k, 1 min, 4°C).
- 24. Resuspend in 1 ml acetone.

**25.** Store tubes at -20°C. Cells are okay for a few days.

\*add one tablet to 10 ml extraction buffer, store at  $-20^{\circ}$ C in 2 ml aliquots

\*\* can also include protease inhibitors at this stage if desired.

## Notes:

1. Cells in the –T control may lose soluble nucleoplasmic proteins even without detergent. It is useful to fix an aliquot of cells directly at step 1 to avoid this problem.

2. Using 100% ethanol for fixation, rather than using methanol/acetone works fine for some GFP tagged strains.

3. To carry out nuclease digestion control, to establish if chromatin integrity is required for the retention of a specific protein, add micrococcal nuclease to 2.5 units/ml and  $CaCl_2$  to 2 mM at step 20 (+Triton). As a control add EGTA to 10mM as well.

## 2. Mounting cells for fluorescence microscopy for visualization of GFP fluorescence

1. Resuspend cells (from step 27) and transfer ca. 0.2 ml to 1.5 ml eppendorf tube. Spin down (4k, 1 min) and take off most acetone - leave about 20-50  $\mu$ l.

2. Resuspend cells by vortexing and mildly sonicate (optional) to break up cell clumps. Spread about 10  $\mu$ l of cell suspension onto a polylysine coated slide (eg BDH 406/0178/00). Allow acetone to dry, and add ca. 10  $\mu$ l mountant. Use agarose mountant for YFP strains.

### Agarose mountant (which seems to preserve YFP fluorescence of Pol1-YFP)

25/5/05

1 g DABCO (antifade\*) 10 ml 0.1 M Tris-acetate pH 8.5

Dissolve, re pH to 8.5 with acetic acid Make up to 100 ml

Add 0.15 g LMT agarose to 25ml and melt in microwave

Add DAPI to 50-125 ng/ml Store in dark (4°C?)

### To use

1. Melt agarose and cool to 37°C

2. Spread cells in acetone in normal way (use ordinary slide – no need to use polylysine coated)

3. Add 10 ul agarose moutant

4. Drop on coverslip – leave 5 min before viewing to allow time to set.

\* doesn't actually seem to make a lot of difference to the rate of fading

## **DAPI/PBS** mountant

50% glycerol 50% PBS 0.005-0.1 μg/ml DAPI (use low conc for CFP filters otherwise can get bleed through)

## 3. Staining extracted cells with antibody

*This is for use with* <u>*Triton-extracted*</u> *cells after methanol and acetone fixation (i.e. from step 7).* 

1. Resuspend Triton-extracted cells and transfer ca. 0.2 ml to 1.5 ml eppendorf tube. Spin down (4k, 2 min) and take off most acetone - leave about 20-50  $\mu$ l.

2. Resuspend cells by vortexing and mildly sonicate to break up cell clumps. Spread about  $10 \mu l$  of cell suspension on to a polylysine coated 13 mm coverslip.

3. Rinse coverslip in PBS, then incubate in PBSBAL for 30 minutes.

4. Remove PBSBAL, add 20  $\mu$ l primary antibody to coverslip. Incubate in a 'wet box' for at least 1h.

5. Wash in PBSBAL (3 x 5mins).

6. Add 20  $\mu$ l secondary antibody (eg Texas Red conjugated). Incubate in dark 'wet box' for at least 1h.

7. Wash in PBSBAL (3 x 5mins).

8. Drain coverslip well. Mount in DAPI/PBS.

(we have stained extracted cells with anti tubulin (TAT1) antibody - gives good staining of mitotic spindles, although interphase microtubules are not preserved).

# PBSBAL

1 X PBS containing:

100 mM lysine hydrochloride, pH 6.90.01% sodium azide1% essentially fatty acid free BSA (Sigma A0281)

# Cover slips

Soak in 0.1 % poly L lysine (5-15 mins). Drain, dry, rinse in water, dry.

Useful stocks

Spermine 1M stock

0.348 g for 1ml of 1M stock

Spermidine 1M stock

0.254g for 1ml of 1M stock

HU

1.388 g for 20 ml 1M stock

Trouble shooting problems

1. DAPI false signal problem with cfp – use dapi at low conc eg 5 ng/ml and integrate to get decent signal

2. when doing hu expts can add hu to zym to arrest

3. essential to get good digestion >95% - use at least 50 ul of zym per digest and prewash in DTT

4. use emm for temp shifts – gives lower bg with mcm4 at least?

5. use emm not yes for temp shifts and do a minimal shift!